





APK Technology

YOUR MOST RELIABLE SUPPLIER OF MEDICAL ACCESSORIES



ECG CABLES CATALOG



ECG CABLES

The reusable and disposable ECG cables designed and manufactured by APK are well recognized by European and American clients. Ever since we had the professional R&D team which has rich experience of working with GE/Philips on board, the quality and stability of our ECG cables can be the top quality products in this industry. Meanwhile, our cables has a very attractive pricing strategy to support our long term partners.



One-Piece ECG cable
3 Leads



Disposable ECG Lead Wire



ECG Trunk Cable



Reusable ECG Lead Wire
5 Leads



New Design Clip



Regular Clip



Snap for Neonatal/Infant



Regular Snap



Nihon Kohden Type



VS Type



Din Type



Euro Type



AA Type



LL Type



ASP Type



Conmed Type (the old type we don't produce now)

Specification

- ◆ Conductor Resistance: $\leq 1\Omega$ (defibrillation circuit resistance not included)
- ◆ Insulation Resistance: wire to wire, wire to shield, $> 20M\Omega$, DC 500V
- ◆ Dielectric Withstanding Voltage:
 - Wire to wire, AC1500V, 50Hz, 0.25mA, No breakdown and flashover
 - Wire to external shell, DC5000V, 50Hz, 0.25mA, No breakdown and flashover

Working Conditions

- ◆ Ambient Temperature: $5^{\circ}\text{C}\sim 40^{\circ}\text{C}$
- ◆ Relative Humidity(RH): 10%~95% (non-condensing)
- ◆ Atmospheric Pressure: 70kPa~106kPa

Cleaning and Disinfection

1. Please unplug the cable from the monitor before cleaning or disinfection
2. Clean or disinfect the cable before attaching to a new patient
3. Clean or disinfecting the cable periodically
4. Clean the cable and patient contact surface with soft cloth moistened with soft dry cloth
5. We recommend using the 75% medical alcohol as the disinfection solution to disinfect for cable and patient contact surface

YOUR MOST RELIABLE SUPPLIER OF MEDICAL ACCESSORIES



APK Technology Co., Ltd.

Shenzhen Headquarter

601, C6 Building, Industry of Hengfeng, Hezhou, Xixiang,

Bao'an District, Shenzhen China

Tel: 86-755-27325581/27325583

Fax: 86-755-27325585

Email: info@apk-technology.com

www.apk-technology.com

